



# WHILL Model F



# The Perfect Travel Power Chair Has Arrived.

Lightweight, foldable, maneuverable. The Model F allows you to explore your world with greater mobility and independence.

Experience the same performance and smart technology that WHILL is known for with a new foldable design.




## Foldable Design

Easily transport the Model F on buses, trains, and planes – or in the trunk of a car or taxi.\*



When not in use, the Model F can be conveniently stored in a closet or out-of-the-way corner of the room.



\*Model F is not to be used as a seat while it is being transported.

The information contained herein is correct at the time of printing, and we reserve the right to alter specifications without prior notice.

## Lightweight

Weighing only 24 kg., the Model F is one of the lightest power chairs on the market today.\*\*



Top Speed



Range




Obstacle Clearance



780 mm turning radius for tight turns

## Connect Intelligently

The WHILL app allows you to check key device information and drive the Model F remotely.



Drive remotely, lock/unlock, and manage speed settings via the WHILL smartphone app.

## Drive with Peace of Mind

Stay engaged in activities and embrace your surroundings knowing the Model F will keep you safe while you ride.



Safely drives up or down 10° slopes.



\*\*Weight is 26.7 kg. if battery is not removed. Battery weighs 2.7 kg.

Note: Driving range of 20 km will vary due to user weight, drive surface, terrain, ambient temperature, speed setting, and battery charge and condition.

# Customization

Choose from five colors and a great selection of accessories to match your personal style and enhance your driving experience.

## Color Options




- White
- Black
- Light Green
- Red
- Light Blue



Controller can be positioned on either the left or right side.

## Accessories



## Optional parts



Dust Cover

# Specifications

Driving Range	Up to 20 km	Maximum Speed***	6 km/h
Foldable	Yes	Max Incline	10°
Turning Radius	780 mm	Obstacle Clearance	35 mm
Seat Sizes Available	18"	Width Between Armrests	Adjustable: 2 positions (17", 19")
Device Width	555 mm / 605 mm	Device Length	935 mm (465 mm when folded)
Device Weight	26,7 kg	Charging Time	~ 5 hours
App Enabled	Yes, with remote control	Air Travel Approved	Yes
Max Weight Capacity	Up to 115 kg	Swappable Arm Covers	5 Colors
Battery	Lithium-Ion (25.3V 10Ah)	Water Resistance	IPX4

\*\*\*Acceleration and/or deceleration may change when driving on steep slopes.

- This device is intended to provide mobility to individuals who are comfortable with operating a power chair.  
- WHILL Model F is classified as Wheelchair Class A.

WHILL Europe B.V.

Johan Cruijff Boulevard 65, 1101 DL, Amsterdam, The Netherlands  
info@whill-europe.com www.whill.inc

